

# iCheck Carotene: Measurement Step-by-step



# Outline

- 1. Introduction to iCheck Carotene**
- 2. Using iCheck Carotene Step by Step**
- 3. Handling and Disposal**



# iCheck Carotene is a Portable All-Inclusive Test Kit



- iCheck Carotene measures carotenoids in vitamin premix, food and biological fluids
- It is a portable single-wavelength photometer, precalibrated for quantitative measurement of carotenoids.
- Carotene reagent vial has a precise volume and ratio of organic solvents for optimal extraction of carotenoids from the sample.




# iCheck Carotene Measures Carotenoids in Food and Biological Fluids

- iCheck Carotene validated to measure carotenoids in following samples:



Beverages



Salmon, eggs



Vitamin Premix



Cattle Blood




Roots (i.e. cassava)



Breast Milk



# Science Behind iCheck Carotene



# iCheck Carotene is Precalibrated and Comes with a 2-year Warranty



- iCheck Carotene is pre-calibrated during production. Therefore, no user calibration is necessary.
- iCheck Carotene comes with a 2-year warranty.
- The calibration of the device is controlled on 2 levels:
  - Automatic Self-Test
  - Device Control with Carotene Standard
- The shelf-life of the reagent vials is 1 year at 20-30 °C.



# Provitamin A carotenoids and conversion to vitamin A

## USA and Canada

1 mg Retinol Activity Equivalent (RAE)

- = 1 mg retinol
- = 2 mg supplemental beta-carotene
- = 12 mg dietary beta-carotene
- = 24 mg other dietary provitamin A carotenoids

## Europe


1 mg Retinol Equivalent (RE)

- = 1 mg retinol
- = 6 mg beta-carotene
- = 12 mg other provitamin A carotenoids

If you need support with conversion, you can use Units Converter available on [iCheckAcademy.org](https://www.icheckacademy.org)



# iCheck Carotene Measurement Range is 0.15 – 15.0 mg/L





- iCheck Carotene results are displayed in mg per litre and indicates the total carotenoids concentration in the sample
- iCheck Carotene measurement range is:
  - 0.15 – 15.0 mg/L

Measurement Limit	Device Display	Observed in the vial
< 0.15 mg/L	Measured Value	No yellow color
> 15 mg/L	Value above 15 mg/L	Intense yellow color






# iCheck Carotene is validated and the results published in peer reviewed journals





- Comparison to the reference method of HPLC yields correlation coefficient of **0.94 to 0.99**.
- Coefficient of variation in cattle blood is **3.5%**.
- *References:*
  - Schweigert FJS. Determination of  $\beta$ -carotene in whole blood of cattle Comparison of a new cow-side assay with HPLC. *Vet. Clin. Path.*, 2011.
  - Islam KMS, Schweigert FJS. Comparison of three spectrophotometric methods for analysis of egg yolk carotenoids. *Food Chemistry* 172, 233–237, 2015.



# iCheck Carotene Measures in a 4-step Procedure



# Step 1: Device Control



## Control your iCheck Carotene:

- Measurement must be done on a flat surface. Keep the metal cap on the measurement chamber.
- Turn on your iCheck™ Carotene. It will run automatic **self-test**.
- Insert Carotene Standard into the measurement chamber.
- Go to **Device Control** measurement mode and press measurement key.
- Control the value displayed against the range indicated on the standard case.
- *Control your device before each set of measurements*


## Attention

- Do not measure outside the temperature range of 20 – 30 °C
- Do not allow any dirt or dust to enter the measurement chamber



# Step 1b: Sample Preparation

For more detailed instructions on specific sample preparation protocols refer to the “iCheck Carotene\_Sample Preparation” file.




- Prepare the sample by diluting it with water.
- Sample preparation should be done if
  1. expected concentration of your sample is above iCheck Carotene measurement range ( $> 15.0$  mg/L)
  2. your sample is a solid sample.

## Important!

- Total Carotenoids concentration of the sample solution must be in the measurement range of iCheck Carotene which is  $0.15 - 15.0$  mg/L.



# Step 2: Injection



## Inject the diluted sample


- Vigorously shake the sample solution to homogenize.
- Quickly take up approx. 0.6 mL of the sample solution with the syringe.
- Put on the needle and make sure no air bubbles are in the syringe.
- Adjust the volume to exactly 0.4 mL into paper tissue.
- Inject 0.4 mL sample into the reagent vial.

## Attention

- The sample must not settle down in the solution
- There should be no air bubbles in the syringe



# Step 3: Extraction



## Extraction of total carotenoids:

- Vigorously shake the vial for 10 seconds. The content of the vial should appear as one uniform solution.
- Let the vial stand still for 5 minutes until the solution in the vial appear as two distinct phases.
- Make sure nothing is stuck on the inside of the glass in the vial on the level of the upper clear phase. If a piece of sample is stuck on the inside tap gently the vial on the solid surface to shake off the piece and wait again for phase separation.


**Note:** some samples do not separate into 2 phases. To facilitate this use BioAnalyt centrifuge.

## Attention

- Do not measure a reagent vial with NO clear supernatant of at least 0.8 cm (see picture)
- Do not store vial in the fridge



# Step 4: Measurement



## Measure the reagent vial with the sample:

- Make sure the vial is clean. Wipe the glass surface with paper tissue
- Set the iCheck Carotene to **OTHER SAMPLE** mode.
- Press measurement key and follow the instructions on the display.
- iCheck™ Carotene will display the result in mg/L. This value indicates the concentration of total carotenoids (TC) in the injected sample solution.
- The vial with the sample must be measured with iCheck earliest 5 minutes and latest 1 hour after sample injection.

## Attention

- Do not measure in the direct sunlight
- Always measure in the temperature range of 20 – 30°C



# Calculations of Dilution Factor and Measured Concentration

1

$$\text{Dilution Factor (DF)} = \frac{\text{Total Diluted Sample Volume [mL]}}{\text{Sample [g]}}$$

2

$$\text{Measured Total Carotenoids [mg/kg]} = \text{iCheck Carotene result [mg/L]} \times \text{DF}$$





# Disposal Instructions for Vials and Syringes




## REAGENT VIALS

Dispose of vials as hazardous material according to local regulations and upon consultation with local waste disposal services. Handle any broken glass using protective gloves.

## SYRINGES

Take extra care when handling syringes with needles. **DO NOT RECAP** the needle!

Discard the needles into specific containers.



# For Technical Support email:

[support@bioanalyt.com](mailto:support@bioanalyt.com)

or WhatsApp:



+49 332 835 150 034

## Contact for orders:

E-mail: [contact@bioanalyt.com](mailto:contact@bioanalyt.com)

[www.bioanalyt.com](http://www.bioanalyt.com)

Follow us on:  